GLOBAL ONLINE AUCTION



Global Online Auction

Late Model Laboratory and R&D Equipment
Formerly of A Palo Alto Pharmaceutical Company

Starting: July 11th, 2016 (7am PDT) | Ending: July 12th, 2016 (10am PDT)



For more information including a photo gallery, equipment list, and more, please visit

www.hgpauction.com

Items subject to prior sale. Buyers are advised & encouraged to inspect equipment and read all Terms and Conditions prior to auction.

LATE MODEL LABORATORY AND R&D EQUIPMENT

EXCESS TO FUTURE NEEDS OF CARBYLAN THERAPEUTICS. BURLINGAME, CA



HPLC SYSTEM (6 MODULES). AGILENT 1100-SERIES (6-MODULES-TOTAL) HPLC SYSTEM. [WITH MODULES: AGILENT 1100/G1365B MULTIPLE WAVELENGTH DETECTOR S/N-DE33202334; AGILENT 1100/G1311A QUATPUMP S/N-DE33224159; AGILENT 1100/G1329A AUTOSAMPLER DE33211350; AGILENT 1100/G1330B AUTOSAMPLER THERMOSTAT S/N-DE33211350; AGILENT 1100/G1316A COLCOM S/N-DE33235135; AGILENT 1100/G1364. ALSO AGILENT 1100-SERIES REAGENT TRAY] [WITH DELL PRECISION-T3400 COMPUTER, 17" MONITOR AND KEYBOARD]



AGILENT 708-DS DISSOLUTION APPARATUS (AGILENT MODEL#G7910A P/N-G7910-60423; S/N-MY15278304) (8-POSITION CONFIGURATION/BUT 6-POSITION TESTING; 1L, MOTORIZED LIFT) (WITH AGILENT #G7986A PUMP). WITH DELL OPTIPLEX-330 COMPUTER, 19" MONITOR AND KEYBOARD



HPLC SYSTEM (7 MODULES): AGILENT 1290/1260/1200/1100 (7-MODULES-TOTAL) HPLC SYSTEM. [WITH MODULES: AGILENT 1200/G1362A REFRACTIVE INDEX DETECTOR S/N-CN60555595; AGILENT 1260/G1365D MULTIPLE WAVELENGTH DETECTOR VL S/N-DEAAZ01006; AGILENT 1100-SERIES/G1311A QUATERNARY PUMP S/N-DE33223865; AGILENT 1100-SERIES/G1322A VACUUM DEGASSER S/N-JP32750041; AGILENT 1100-SERIES/G1329A AUTOSAMPLER S/N-DE91605034; AGILENT 1290-INFINITY/G1330B AUTOSAMPLER THERMOSTAT S/N-DEBAK21905; AGILENT 1100-SERIES/G1316A COLUMN COMPARTMENT S/N-DE23931977. ALSO AGILENT 1200-SERIES REAGENT TRAY. INCLUDES AGILENT 1200-INFINITY SERIES TOOLS KIT] [WITH DELL OPTIPLEX-360 COMPUTER; DELL 24" MONITOR; KEYBOARD AND MOUSE]



PLC SYSTEM (6 MODULES): AGILENT/HP 1100-SERIES (6-MODULES-TOTAL) HPLC SYSTEM. [WITH MODULES: AGILENT 1100/G1315A DIODE ARRAY DETECTOR S/N- DE91608025; AGILENT 1100/G1312A BINARY PUMP S/N-DE91605415; AGILENT 1100/G1329A AUTOSAMPLER S/N-DE11106926; AGILENT 1100/G1330A AUTOSAMPLER THERMOSTAT S/N-DE82207677AGILENT 1100/G1316A COLCOM S/N-DE03016832; HP 1100/G1322A VACUUM DEGASSER S/N-JP73021018. ALSO AGILENT 1100-SERIES REAGENTTRAY] [WITH DELL OPTIPLEX-GX150 COMPUTER, 15" MONITOR AND KEYBOARD].





GLOBAL ONLINE AUCTION: JULY 11TH, 7AM - JULY 12TH, 10AM PDT PREVIEW: JULY 11TH 9AM - 4PM OR BY APPOINTMENT - CONTACT AUCTIONEER FOR FURTHER INFORMATION



HPLC SYSTEM (6 MODULES): AGILENT 1290/1260/1100 (6-MODULES-TOTAL) HPLC SYSTEM. [WITH MODULES: AGILENT 1260/G1315B DIODE ARRAY DETECTOR VL S/N-DEAAX06706; AGILENT 1260-INFINITY/G1312B INFINITY BINARY PUMP S/N-DEACB08937; AGILENT 1260/G1329B AUTOSAMPLER S/N-DEAAC30798; AGILENT 1290-INFINITY/ G1330B AUTOSAMPLER THERMOSTAT S/N-DEBAK21904 ; AGILENT 1100-SERIS/G1322A VACUUM DEGASSER S/N-JP7302310 ; AGILENT 1100-SERIES/G1316A COLUMN COMPARTMENT S/N-DE91609599. ALSO AGILENT 1200-SERIES REAGENT TRAY. INCLUDES AGILENT 1200-INFINITY SERIES TOOLS KIT] [WITH DELL OPTIPLEX-755 COMPUTER; DELL 19" MONITOR; KEYBOARD AND MOUSE]



MILLROCK TECH LD85-B3 LYOPHILIZER S/N-LD850906007B3. WITH LEYBOLD TRIVAC **#D8B VACUUM PUMP**



LOT: THALESNANO ICECUBE FLOW REACTOR SYSTEM. NEVER USED. STILL IN BOX. CONTENTS OF (4) BOXES TO INCLUDE: ICECUBE REACTOR MODULE, PUMP MODULE, AND CONTROL UNIT (PRODUCT CODE THS 09030). THE REACTOR HAS 2 REACTION ZONES WHICH CAN BE CONNECTED FOR MULTISTEP REACTIONS. THE FIRST REACTION ZONE HAS A TEMPERATURE RANGE -70° C OR - 50°C (SPECIFIED BY USER ON PO) AND +80°C. THE FIRST REACTION ZONE HAS A MAX VOLUME OF 8ML (-50°C VERSION) OR 4ML (-70°C VERSION). THE SECOND REACTION ZONE HAS A TEMPERATURE RANGE OF -30 - +80°C, AND SUPPORTS A REACTION VOLUME OF 4ML, THE REACTOR MODULE IS USED WITH 1/16" OR 1/8" TEFLON TUBES AND CAN BE PRESSURIZED UP TO 10 BAR. THE PUMP MODULE HAS TWO BUILT IN ROTARY PISTON PUMPS. EACH PUMP HAS FLOW RATES OF 0.2 - 0.4 ML/MIN. BOTH PUMPS ARE CONNECTED VIA 1/8" PTFE TUBES TO THE INLET SWITCH VALVES THAT CONTROL EACH OF THE INPUT LIQUID SOURCES. THE CONTROL UNIT IS A USER INTERFACE THAT CONTROLS ALL PARAMETERS OF A REACTION, AND ALL ASPECTS OF THE CONNECTED MODULES WHICH ARE PART OF THE ICECUBE SYSTEM. THIS UNIT INCLUDES (2) ICECUBE SINGLE PUMP MODULES (PRODUCT CODE THS 09032) THAT FUNCTION WHEN MORE THAN 2 LIQUID FEEDS ARE REQUIRED FOR A REACTION. THIS SEPARATE UNIT CAN BE CONTROLLED INDEPENDENTLY USING FLOW RATE RANGE OF 0.2ML - 0.9ML/ MIN ON THE TOUCH SCREEN OF THE PUMP MODULE OR VIA THE CONTROL UNIT.



SP SCIENTIFIC MODEL-VIRTIS 35L GENESIS SQ EL-85, LYOPHILIZER WITH ENCORE CONTROL SYSTEM (1200L, 6 SHELF BULK, CLEANROOM NARROW CABINET STYLE). ITEM STILL 'NEW IN BOX' AS WRAPPED/PACKAGED FOR SHIPPING BY MANUFACTURER. [LYOPHILIZER: P/N-100010546; S/N-320549; YEAR-2016; 208-230 VOLTS 1-PHASE 30-AMPS 60-HZ] [INCLUDED: GENESIS MULTI-LEVEL BILL OF MATERIALS BUILD-UP STANDARD REPORT; ENCORE MULTI-LEVEL BILL OF MATERIALS BUILD-UP STANDARD REPORT; DRAWINGS/SCHEMATICS; VALIDATION GUIDE AND WORKBOOK MANUAL: CALIBRATION AND CONFORMITY CERTIFICATES/DOCUMENTATION: QC TESTING DATA REPORT; LYOPHILIZER OPERATIONS MANUAL] [PACKAGED DIMENSIONS APPROX: MAIN UNIT 76"LX35"WX81"H ON CASTERS; BOX#1 25"LX25"WX38"H; BOX#2 15"LX14"WX9"H]





AUCTIONEER'S NOTE

HGP is pleased to be selected to conduct this Global Online Auction for:

Carbylan Therapeutics

ABOUT THE SALE

This sale features Late Model Laboratory and R&D Equipment Formerly of A Palo Alto Pharmaceutical Company.

Please view our website at **www.hgpauction.com** for the bidding instructions and the full terms of the sale. Thank you for your interest, and we look forward to providing you quality equipment sales and World Class Service.

KEY ASSETS

Thalesnano Flow Reactor • Freeze Dryer • Millrock LD85 Lyophilizer • Agilent 1100/1200 Series HPLC Systems • Agilent 1260 1362A RID • Stereoscopes • Microscopes • Autoclaves Rheometers • Environmental Specialty Chambers • Viscometer • Waterbaths • Balances Shakers, Stirrers & Hot Plates • Incubators and Oven • Glassware Consumables • And More



Scan this QR Code with your smartphone for More Information.

Contact Information

Nick Dove, Tel: (858) 847.0659 Email: ndove@hgpauction.com

A BUYER'S PREMIUM WILL APPLY AT THIS SALE

For more information including a photo gallery, equipment list, and more, please visit us online at www.hgpauction.com. Items subject to prior sale. Buyers are advised & encouraged to inspect equipment and read all Terms and Conditions prior to auction.



GLOBAL ONLINE AUCTION: 400 AUCTION LOTS

Late Model Laboratory and R&D Equipment Formerly of A Palo Alto Pharmaceutical Company

GLOBAL ONLINE AUCTION

Starting: July 11, 2016 – 7:00 am PDT **Ending:** July 12, 2016 – 10:00 am PDT

AUCTION LOCATION

1576 Rollins Road, Burlingame, CA

AUCTION PREVIEW

Monday, July 11, 2016 (9:00 am - 4:00 pm PDT) Contact Nick Dove (858) 847.0659 to schedule a private preview.

For more information including a photo gallery, equipment list, and more, please visit us online at www.hgpauction.com

Items subject to prior sale. Buyers are advised & encouraged to inspect equipment and read all Terms and Conditions prior to auction.

WOULD YOU LIKE TO SELL AN ASSET?

Simply call 1.650.207.0627 and speak with Kirk Dove, Managing Partner for more information. We will be glad to listen to your specific needs and from that HGP can prepare a detailed solution for your review that will maximize the returns for you're assets for the company and It's Shareholders.

GLOBAL HEADQUARTERS

Hacienda Del Mar 12625 High Bluff Drive San Diego, CA 92130 Phone: (877) 303.8040 Fax: (866) 388.8595

FIRST CLASS MAIL

POSTMASTER: Dated Material. Please Deliver Promptly.

Presorted First Class Mail U.S. Postage PAID San Diago, CA Permit No. 00

Starting: July 11, 2016 – 7:00 am PDT Ending: July 12, 2016 – 10:00 am PDT Location: 1576 Rollins Road, Burlingame, CA

For complete Terms of Sale, visit our website at:

www.hgpauction.com



















